## JUP ITER SERIES 500

AUTOMATIC GAS CHROMATOGRAPHIC ANALYZER for BTX and  $\mu$ VOCs.

NIRA Analyzer Mod. **JUPITER 501** is a very sensitive instrument specifically designed to automatically measure traces of organic volatile components in air.

It is an automatic equivalent system for sampling and measuring:benzene, toluene, ethyl-benzene and xylene [Analytic method of reference according to DM on 25/11/1994] or Volatiles Organic Composites.

It was designed to simplify typical monitoring problems, carrying out the task in fixed cabins or mobile units, maintaining the sensitivity and selectivity of typical gas chromatographs in laboratories.

This instrument determines organic volatile composites "enriching" a trap containing absorbent material, then "desorbing" it in a capillary column for detection via FID.

With the new **JUPITER 501**, N.I.R.A. gives coherent and complete results to those who value air quality: urban vehicular traffic pollution monitoring in fixed or on mobile units; industrial or residential area monitoring near industrial plants that may produce aromatic components.

A specific Jupiter model system is suitable to monitor the air quality in clean rooms.



Jupiter 500

## BENEFITS

- Flame Ionization Detector: enhanced linearity, reproducibility, reliability and can be used in cabins.
- > Specific sampling and switching valves: they can work for years at a high temperature without any maintenance.
- > Optimized pneumatics: an automatic electronic pressures and temperatures monitoring system is incorporated in a traditional pneumatic circuit, this allow easier interface and performances optimization.
- Sampling using a pump and a selectable mass-flow controller, all integrated in the system.
- > Enrichment on trap: it does not require added cryofocusing systems.
- Automatic gas chromatographic separation with long life capillary columns.
- Automatic injection system and thermal desorption.
- Hermetic detectors with channeled vent: they can be used in any area, even in controlled atmospheres.



## "PC EMBEDDED" TECHNOLOGY

- Allows the total visibility in every analytical step.
- Graphic and numeric data re-loading off-line.
- Increase the quality of chromatograms thanks to an LCD/TFT display and a real time analysis control.
- Continuous checking of: operating parameters, analysis and conditioning.

• Trends recordings.

- Increased sampling accuracy and precise chromatogram elaboration.
  - te data elaboration: integration, trend Auto-diagnostic.
- Accurate data elaboration: integration, trend recording, and history of all anomalies that may occur.

## ■ TECHNICAL SPECIFICATION

ANALYZED COMPONENTS	Benzene, toluene, ethyl-benzene, xylene	
RANGES	0-1000 ppb F.S.	
RESPONSE TIME	15 minutes (30-60 minutes = optional)	1
NOISE	0,1 ppb	
MINIMUM DETECTABLE	0,3 ppb	
SAMPLING FLOW AND TIME	Selectable	
MAXIMUM SAMPLING VOLUME	1000 ml	
ACCURACY	1%fs	
LINEARITY	1%fs	
OPERATIONAL TEMPERATURE	+5; +40°C	- 5
ANALOGUE OUTLETS	(0-10)V, (4-20)mA	
MAXIMUM NUMBER OF ANALOGUE OUTPUTS	8	
SERIAL OUTPUT	RS232	- 5
PRINTER OUTPUT	USB	
DISPLAY	LCD - 10.4"	- 1
KEYBOARD	PC compatible	
PROGRAMMABLE SAMPLE INTAKE CAPACITY	50-300 ml/min.	i i
CALIBRATION	Local and remote	- 1
HYDROGEN CONSUMPTION	60 ml/min.	
AIR CONSUMPTION	300 ml/min.	
HYDROGEN CONSUMPTION	2 bar minimum	- 1
AIR PRESSURE	4.5 bar minimum	
DIMENSIONS	450x270x600mm (L x W x H)	
WEIGHT	30 kg	
POWER	115/230 Vca, 50/60 Hz, 300W	- 1
TO STATE OF THE PARTY OF THE PA	The state of the s	100



140