

# 5380 PFPD

## PULSED FLAME PHOTOMETRIC GC DETECTOR

OI Analytical's patented\* 5380 Pulsed Flame Photometric Detector (PFPD) excels at selective, high sensitivity detection of sulfur, phosphorus, and other elements. The 5380 PFPD provides a 10-100x increase in signal-to-noise and a 10x increase in selectivity over traditional Flame Photometric Detectors and requires significantly less maintenance and gas to operate.

The 5380 PFPD consists of a 5380 Detector assembly, 5380 Detector Controller, pneumatic components located either in the Detector Controller or in the gas chromatograph (GC), and WinPulse user interface software.



### Operating Principle

The PFPD operates using a propagating flame that terminates within a quartz combustor. Gas phase reactions produced by the propagating flame result in light emissions with specific luminescent spectra and lifetimes. The differences in specific emission lifetimes combined with the kinetics of the propagating flame allow both time and wavelength information to be used to improve the PFPD's selectivity and to decrease the observed noise, enhancing sensitivity. The propagating flame uses low combustible gas flow rates, increasing the relative analyte concentration. In addition, use of gated electronics permits acquisition of two simultaneous, mutually selective chromatograms.

### PFPD Capabilities

- Superior sensitivity and increased selectivity for S and P compared to conventional FPDs
- Linear, equimolar response for quick, easy calibrations
- Simultaneous mutually selective chromatograms e.g.: S+C, or S+P
- Inherent self-cleaning design completely eliminates soot formation, or "coking", seen in other sulfur-selective detectors
- Long-term stability and minimal maintenance

### Principal Applications

- Sulfur in petrochemicals
- Organophosphorus pesticides
- Flavor and fragrance analyses
- Sulfur in beverage-grade CO<sub>2</sub>
- Simultaneous PFPD and MS detection
- Chemical warfare agents
- Organotin compounds in environmental samples
- Organometallic detection
- Explosives analysis
- P, S, As, and Si detection in the semiconductor industry
- Sulfur in pharmaceuticals

# 5380 PFPD SPECIFICATIONS

## Performance Specifications

<b>Detectivity</b>	
Sulfur	<1 pg S/sec
Phosphorus	<100 fg P/sec
<b>Sensitivity</b>	
Sulfur Signal-to-Noise	>300 (at 10 pg S/sec elution rate peak-to-peak noise)
<b>Drift (S or P)</b>	<10x peak-to-peak noise in 20 min
<b>Selectivity (at Optimum Detectivity Levels)</b>	
Sulfur	> 10 <sup>6</sup> S/C
Phosphorus	> 10 <sup>5</sup> P/C (selectivity is adjustable with a trade-off in detectivity)
<b>Detector Linearity</b>	
Sulfur	Quadratic in response. Linear to approximately three orders of magnitude.
Phosphorus	First order linear over approximately three orders of magnitude.
<b>Response Uniformity</b>	Equimolar ± 8% (S, P)
<b>Chromatographic Peak Tailing</b>	<0.2 sec in S and P
<b>Gas Requirements</b>	
Carrier	He or H <sub>2</sub> at 80 psig; 99.8% purity or better
Air	60 psig; zero air (CGA grade E)
Hydrogen	60 psig; 99.995% purity or better (electrolytic grade)
<b>Power Requirements</b>	115/230 V <sub>AC</sub>
<b>Computer Requirements</b>	
Operating System	Windows® XP Pro, Windows® 7
Comm Ports	One serial (RS-232), 16550 UART or USB to Serial Adaptor
<b>Minimum Temperature</b>	180 °C
<b>Maximum Temperature</b>	420 °C
<b>Carrier Gas</b>	5 mL/min maximum flow rate Helium; up to 10mL/min using H <sub>2</sub> carrier gas
<b>Typical Gas Consumption</b>	
H <sub>2</sub>	10-15 mL/min
Air	20-30 mL/min
<b>Humidity</b>	5-80% relative humidity
<b>Altitude</b>	2,000 m maximum

<b>Controller Board Inputs and Outputs</b>	
Two Channels	0-1 V
One Serial	RS-232-C
One Signal In	Electrometer; PFPD
High Voltage Out	PMT 0-1,000 V
Ignitor Current	0-3.4 A
Oscilloscope Output	20 Hz, 25 ms display
S/W HV Protection	PMT Protection
Timed Events (from GC Remote Start)	Autozero, range, attenuation, ignitor, mode or channel (e.g. S, P, C), and record
<b>Controller Dimensions</b>	22.2 cm H x 14 cm W x 33 cm D (8.75" H x 5.5" W x 13"D)
<b>Pneumatic Control (Standard)</b>	Manual control of detector gases with mass flow controllers and metering valve

## Options

<b>PFPDView Software</b>	For post-acquisition signal processing
<b>Pneumatic Control</b>	Automatic using GC electronic flow control of detector gases
<b>Combustor</b>	3 mm

## Simultaneous sulfur and hydrocarbon chromatograms of gasoline containing 5 ppm total sulfur.

